

Product Name :
Universal Testing Machine Rotating Fatigue Machine

Product Code :
ASTRALAB-SUPPLYA210100



Description :

Universal Testing Machine Rotating Fatigue Machine

Technical Specification :

Autographic recorder: Simultaneous roll autographic recorder supplied as standard to enable study of the behaviors of materials.

Ideal Dia.: High reading accuracy due to large size ideal design of dial.

Large columns: Large effective clearance between columns enables testing of standard specimen as well as structures.

Facilities for tests: Motor-driven threaded columns for quick and convenient adjustment of lower cross head to facilitate rapid fixing of tear specimen.

Easy Changeability: Easy change from plain to threaded and screwed specimens.

Simple and Safe: Simple to operate. Robust construction. Chrome plated metal components.

Wide range of standard and special attachments / accessories available.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265