


Product Name :
Nurse Call Trainer

Product Code :
EL - CT3



Description :

Nurse Call Trainer is a trainer on a wall-mounted panel. It is built so that students can assemble, wire, and connect different circuits. The visual nurse call panel manages calls initiated at patient stations and staff duty stations on a central panel that provides visual and tone annunciation of call events. Edwards components are used for this trainer. The trainer comes with a complete reference manual.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265