

Product Name :
Aerodynamics Laboratory Equipment Laboratory

Product Code :
ASTRALAB-SUPPLYA56001



Description :

Aerodynamics Laboratory Equipment Laboratory

Technical Specification :

Schlieren Apparatus For Supersonic Wind Tunnel

Pitot Static Traverse

Bench Top Wind Tunnel

Modular Air Flow Bench

Variable Speed Two Stage Compressor Study Unit With Digital Torque Indicator

Nozzle Flow Apparatus

Subsonic Wind Tunnel Starter Set

Set Of Sectioned Components For Gas Compression

Centrifugal Fan Study Unit

Basic Lift And Drag Balance

Compressors Test Bench

Differential Pressure Unit

32 Way Pressure Display Unit

Centrifugal Compressor Module

Subsonic Wind Tunnel

Double Stage Compressor Study Unit

Axial Fan Module

Variable Scale Differential Micromanometer

Fluid Threads Visualizer

Aerofoil With Tappings

Wind Tunnel Flight Demonstration

Compressors Study Unit Reciprocating And Rotative

Set Of Sectioned Components For Pneumatic Controls

Flow Visualization Wind Tunnel

Centrifugal Compressors Study Unit



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265