

Product Name :
OSBORNE REYNOLDS APPARATUS

Product Code :
H25



Description :

Horizontal test pipe diameter 14 mm with bell mouth, pipe length 720 mm. , transparent test pipe shell with stilling baffle, 0.5 lt dye reservoir with a needle. 2 lt measuring cup transparent tank.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265