

**Product Name :**  
Network Analyzer

**Product Code :**  
ASTRALAB-SUPPLYA132009



**Description :**

Network Analyzer

**Technical Specification :**

Two-port vector network analyzer with a full S-parameter test set for bidirectional measurements

Wide dynamic range of typical 120dB

Measurement bandwidths from 1Hz to 500KHz

Wide 10.1" WXGA multitouch screen

System keys for setup, presets, settings

Clearly structured user interface

Undo/redo softkey for user-friendly operation

Fully integrated context-sensitive help menu

Calibration units for quick calibration

De/embedding functionality and fixture compensation



## Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265