

Product Name :
RF Signer Generator Trainer

Product Code :
ASTRALAB-SUPPLYA132010



Description :

RF Signer Generator Trainer

Technical Specification :

Wide frequency range: 10-600MHz.

Internal/external modulation (AM/FM) Fine control over Frequency /
Amplitude range Bright 4 1/2 digit LED/16X2 LCD display for freq.

Two models- manually operated or micro-controller based.

Useful for testing and measurement in RF laboratory.

Optionally PC GUI based control using USB IO Module (UCOM).



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265