

**Product Name :**  
PH Meter

**Product Code :**  
ASTRALAB-SUPPLYA139004



**Description :**

PH Meter

**Technical Specification :**

Fast and Accurate: No need to work with color charts or messy test strips, this meter get the result fast and accurate. It with 0.01 resolution, 0.05 accuracy, full pH 0-14 measurements! And this ph meter is factory calibrated. And when the product with problems, users can make calibration to get problem solved.

Wide Use: This PH meter is great for a wide range.

Easy calibration: No need any adjustment screwdrivers, just press a button to finish the calibration.

ATC Function: ATC accounts for changes in temperature automatically, it can make the meter more accurate.



## **Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265