

**Product Name :**  
Laboratory Constant-Temperature Drying Oven

**Product Code :**  
ASTRALAB-SUPPLYA141002



**Description :**

Laboratory Constant-Temperature Drying Oven Supplier in Africa and Lab Equipment South Africa, Africa Kenya, Uganda, Sudan, Algeria, Morocco, Angola, Ghana, Mozambique, Madagascar, Ivory Coast, Cameroon, Niger, Burkina Faso, Mali, Malawi, Nigeria, Ethiopia, Egypt, DR Congo, Tanzania, Togo, Sierra Leone, Libya, Congo, Central African Republic, Liberia, Mauritania, Eritrea, Gambia, Botswana, Namibia, Gabon, Lesotho, Guinea Bissau, Equatorial Guinea, Mauritius, Eswatini, Djibouti, Reunion, Comoros, Cape Verde, Western Sahara, Mayotte, Sao Tome, Seychelles, Zambia, Chad, Somalia, Senegal, Zimbabwe, Guinea, Rwanda, Benin, Burundi, Tunisia, South Sudan.

**Technical Specification :**

Laboratory Constant-Temperature Drying Oven 70L capacity

External material: cold-rolled steel with anti-bacteria powder coating.

Stainless steel inner chamber, round angle structure, adjustable shelves.

Double-layer glass observation window, adjustable airtight buckle lock.

PID control with LED display

Over-temperature protection.



## Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265