


**Product Name :**  
MINI FLOW CHANNEL

**Product Code :**  
H28



**Description :**

Channel : 620 mm long x 15 mm wide x 180 mm deep and tank with stilling materials. Models : small cylinder, large cylinder, narrow crested weir, broad crested weir, crump weir, symmetrical aerofoil, asymmetrical aerofoil. 0.5 lt dye reservoir with 5 stream lines.



## **Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265