


Product Name :
Bending Moment

Product Code :
SEENGG01AS



Description :

Bending Moment

Technical Specification :

Bending Moment - To work with the universal test frame. Two-part beam, simply supported, with a ball bearing hinge at the joint. Underhung spring balance with adjustable length holds the hinge in compression, to resist bending moment. Movable load hangers fit around the beam, set of weights.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265