


**Product Name :**  
Breadboard Series 2

**Product Code :**  
- BRD-LABAST-002



**Description :**

Breadboard Series 2, We are manufacturer, supplier, exporter and solution provider in quality Breadboard Series 2

**Technical Specification :**

Breadboard Series 2 If you are interested in this product, please download the PDF for more details.- With considerable system expansion available for the computer (Binding post and various optional connectors are available for end users to fit on connector socket).

- With EP-3 edge plate, it accepts big size electronic part with big size or without pin such as potentiometer, switch, light etc....
- 4 binding posts are equipped for external power supply input.
- Special suggested sturdy jump wire is optionally available for connecting pin of universal connector and breadboard.
- Specially suggested KS-350 jump wire kit is optional and suitable for layout of circuit.

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265



**Astra Scientific**