

Product Name :
Multi-Board (MB-85 & PLCC)

Product Code :
BRD-LABAST-006



Description :

Multi-Board (MB-85 & PLCC), We are manufacturer, supplier, exporter and solution provider in quality Multi-Board (MB-85 & PLCC)

Technical Specification :

Multi-Board (MB-85 & PLCC)K&Hs PLCC test boards offer users a fast, safe and solderless method of prototyping circuits with surface mount chip carriers, eliminate time-consuming soldering and expenses of building permanent connections. The PLCC test boards are composed of a chip carrier socket and special designed pins with high quality PC board. The PLCC test board comes in 5 sizes : 20, 28, 44, 68 and 84 pins.

Combine PLCC test boards with MB-85 multiboard, a convenient and useful tool will be created. The MB-85 multiboard's center channel distance can be adjusted according to IC size. Aside from PLCC test boards, LSI, MSI and DIP size chips can all be built easily.

MB-85

- Easy and fast prototyping with surface mount chip carriers
- Eliminates time-consuming soldering
- Compatible with K&H solderless breadboard

SPECIFICATIONS

- 5 interconnected clip.....128
- 25 interconnected clip.....8
- IC capacity (14 Pin).....9

- Weight (gm).....140
PLCC
Model
A
B
C
D
Weight
PL-20
25.4(1.0")
22.86(0.9")
31.88(1.3")
27.06(1.1")
5.6g
PL-28
27.94(1.1")
33.02(1.3")
34.42(1.4")
37.22(1.5")
7g
PL-44
33.02(1.3")
53.34(2.1")
39.5(1.6")
57.54(2.3")
12g
PL-68
40.64(1.6")
83.82(3.3")
47.12(1.9")
88.02(3.5")
23g
PL-84
45.72(1.8")
104.14(4.1")
52.2(2.1")
108.34(4.3")
30g



Astra Scientific

Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265