

Product Name :
PSB-01 AC/DC Power Supply

Product Code :
BRD-LABAST-007



Description :

PSB-01 AC/DC Power Supply, We are manufacturer, supplier, exporter and solution provider in quality PSB-01 AC/DC Power Supply

Technical Specification :

PSB-01 AC/DC Power Supply WITH SHORT-CIRCUIT PROTECTION REGULATED POWER SUPPLY

- Power switch with lamp
- Just plug in and start to use
- Compact and light-weighted
- For both digital/analog circuits

SPECIFICATIONS

- DC power supply:
 - * DC + 5V, 1 A ; -5V, 500mA. fixed
 - * DC 0 to + 15V, 500mA variable
 - * DC 0 to - 15V, 500mA variable
- AC power supply:
 - * AC voltage: 12V / 6V / 0V / 6V / 12V, 300mA fixed
- Size: 185 x 155 x 66mm(L x W x H) with stand
- Weight: 3 Kg

ACCESSORIES

- AC cable x 1
- DC output wires: Red x 4, Black x1



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265