

Product Name :
Logic Probe

Product Code :
BRD-LABAST-010



Description :

Logic Probe for Electrical Training Labs, We are manufacturer, supplier, exporter and solution provider in quality Logic Probe for Electrical Training Labs

Technical Specification :

Logic Probe - Logic Testing TTL and CMOS, Displaying pulse presence and logic states, Catching pulse to 30ns or pulse trains to 17MHz, Input overvoltage protection, ABS case, Two audio tones, Two external terminals

SPECIFICATIONS

- Input impedance. 1M Ω
- Width of min. pulse: 30 ns.
- Max.input frequency: 17 MHz.
- Measurable logic gate
- * Logic1(red LED lights) and buzzer:
TTL: 2.3V \pm 0.2V
CMOS: 70%Vcc.
- * Logic 0 (green LED lights) and buzzer:
TTL: 0.8V \pm 0.2V
CMOS: 30%Vcc.
- Input overvoltage protection: \pm 220V DC/AC 15sec.
- Pulse(yellow LED) flashing time: 500ms.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265