

Product Name :
Bench grinders

Product Code :
AS3E25



Description :

Bench grinders Grinder 12" 180 Grit, 3600 rpm;
- Two 12" Platen (120 or 180 grit,
variable speed

The compact and Bench Grinder Machine for stationary grinding, deburring and finishing operation.
Absolutely different with a standard bench grinder. Ideal for 200mm Out dia and width up to 50mm valgro
deburring and finishing wheels.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265