

Product Name :
Furnace

Product Code :
ASA666



Description :

Furnace For smelting both ferrous and non ferrous metals. Oil or coal fired 50kg molten metal capacity Universal precise laboratory electric furnace designed for materials heat treatment, ceramic and stoneware samples firing. The furnaces could be equipped with folding or side open door. Furnace can be used in laboratories, educational institutions, ceramic studios, and industrial laboratories.

- T max. 1100°C
- Vacuumed ceramic fibre chamber
- Heating elements are integrated around a chamber
- High quality thermoinsulation materials
- Micro processor temperature controller
- Excess temperature controller
- Door interlock switch
- Ceramic bottom plate
- Short heating time
- Hermetic closing door
- Low power consumption
- High precision level



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265