

Product Name :
Creep Test Machine

Product Code :
SM1006



Description :

Creep Test Machine Creep Machine, SIZE: 570 mm x 220 mm,
Operating temperature range: +5°C to +40°C,
Sound Levels: Less than 70 dB,
SPECIMENS: 10 x CP1010 Lead, 10 x CP1020 Polypropylene.
Test weights : 3 x 500 g, 2 x 200 g, 1 x 100
versatile Data Acquisition System (VDAS-B),
Connectivity Kit (SM1006CK) for connection to VDAS®



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265