

Product Name :
Fatigue Test Machine

Product Code :
WP 1401



Description :

Fatigue Test Machine Fatigue Testing Machine, Motor: speed: 2800rpm, output: 0.37Kw.

Load: 0 -300N, Load cycle counter: electronic,
8-digit digital display, can be switched to display speed,

Test bars: tempering steel Ck 35, 230V , ~50Hz



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265