

**Product Name :**  
TILTING FLOW CHANNEL, 100 mm. wide

**Product Code :**  
FM156



**Description :**

10 cm wide x 30 cm high x 400 cm long &ndash;2 to +5% tilting, with pump, feed tank, storage tank accessories and models. Built in adjustable bed at mid section, adjustment range  $\pm 10$  cm



**Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265