

**Product Name :**  
Centrifugal Force Apparatus

**Product Code :**  
TM 600



**Description :**

**Centrifugal Force Apparatus Specification**

- [1] Measurement of centrifugal force on rotating masses
- [2] 5 different path radii
- [3] 3 different rotating masses
- [4] Force measuring range 0...25N
- [5] Rotational speed continuously regulated
- [6] Digital display of force and rotational speed
- [7] l x w x h 420x400x270mm

**Technical Data**

Path radii: 25mm, 50mm, 75mm, 100mm, 125mm

Rotating masses: 54g, 79g, 105g

Speed range: 0...500rpm

**Load**

- measuring range: 0...25N

- resolution: 0.1N

**Dimensions and Weight**

l x w x h : 420 x 400 x 270 mm

Weight : approx. 23 kg

**Connections**

230V, 50Hz, 1 phase

**Scope of Delivery**

1 complete measuring unit



**Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265