

**Product Name :**  
Inertia in Rotational Motion Apparatus

**Product Code :**  
TM 610



**Description :**

Inertia in Rotational Motion Apparatus Rotating rod: 550mm long, D=10mm  
Weights: 100g, 200g, 400g  
Solid cylinder  
- diameter: D=120mm  
- mass: 900g  
Hollow cylinder  
- outer diameter: D=120mm  
- inner diameter: d=110mm  
- mass: 0.9kg  
Weight for the drive: 1N



## Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265