

Product Name :
TILTING FLOW CHANNEL, 150 mm. wide, Adjustable bed.

Product Code :
FM158



Description :

15 cm wide x 35 cm high x 600 cm long, and –1 to +3% tilting channel with pump, feed tank, storage tank, accessories, and models. Built in adjustable bed at mid section, adjustment range ± 10 cm

Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265

