

Product Name :
Moment Of Inertia setup

Product Code :
TM 610



Description :

Moment Of Inertia setup
Inertia in Rotational Motion Apparatus
Technical Data
Rotating rod: 550mm long, D=10mm
Weights: 100g, 200g, 400g
Solid cylinder
- diameter: D=120mm
- mass: 900g
Hollow cylinder
- outer diameter: D=120mm
- inner diameter: d=110mm
Weight for the drive: 1N



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265