

Product Name :
Bell Crank Lever Apparatus.

Product Code :
AA15



Description :

Bell Crank Lever Apparatus. Lever horizontal arm 75cm, vertical arm 20cm, spring balance 10kg capacity, adjustable; slotted cast iron weights, 0.5kg. With one hanger of 0.5 kg and steel mounting.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265