

Product Name :
Boyle's Law Demonstration / PC

Product Code :
WL 102



Description :

Boyle's Law Demonstration / PC COMPRESSOR / VACUUM PUMP: power output: 60W, delivery side: 2bar, intake side: 213mbar, TEMPERATURE MEASURING RANGE: vessel 1: 0 - 80°C, vessel 2: 0 - 80°C, PRESSURE MEASURING RANGE : vessel 1: 0 - 4bar absolute, vessel 2: 0 - 2bar absolute, VOLUME MEASURING RANGE: vessel 1: 0 - 3L. 1 LabVIEW software CD + USB cable



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265