

**Product Name :**  
Condensation Unit

**Product Code :**  
WL230



**Description :**

Condensation Unit HEATER: power: 3kW, adjustable, CONDENSER PIPES: 1x copper, sand blasted, 1x copper, gold plated, surface area: 36.18cm<sup>2</sup>. WATER JET PUMP: required water pressure: 0.5bar, flow rate: 4...12L/min, final pressure: 16mbar. PROCESS TANK: pressure: 2bar absolute. MEASURING RANGES: pressure: 0 - 4bar absolute, flow rate: 2x 0,2 - 6L/min. 230V, 50/60Hz, 1 phase, 1 LabVIEW software CD + USB cable.



## Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265