

Product Name :
Diesel Phasing & Calibration Machine

Product Code :
12PSDW



Description :

Diesel Phasing & Calibration Machine Common/Standard Type Test Bench,

Adjustable rotation speed range: 0~4000 RPM,

Test cylinders: 12,

Rotation speed presetting: 10 kinds,

DC. Supply: 12/24V,

Feeding pump motor: 1.1kw,

Series of graduate glass: 45cc, 150cc. The 260cc,

Volume of oil tank: 60L, 3-phase power supply: 380v/50hz.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265