

Product Name :
Common rail injector tester

Product Code :
CRT-10011



Description :

- Common rail injector tester 1.High pressure 0-1800bar
2. Check the injector starting pressure
 - 3.Assessment of the injector spray fog quality
 4. Inspeicon of the injector spray oil range
 - 5.Check of the injector seal and angel
 - 6.Check of the injector oil flow
 - 7.The machine has fog dissipation convulsion device to improve the work condition.
 - 8.It has the two garde filter to ensure the cleanliness of test.
 - 9.The parameter of common rail injector (fuel spray counts, frquency, pulse width, pressure etc..) can be set willfully.
 10. it also can clean the common rail injector by Ultrasonic Cleaning Instrument.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265