

Product Name :
Digital to Analogue Converter

Product Code :
AEE 9



Description :

Digital to Analogue Converter Four isolated 16-bit outputs

- 100-kHz/channel update rate
- 480 Ksamples/channel max buffer
- ± 1 , ± 2 , ± 5 , and ± 10 VFS programmable unipolar and bipolar output ranges



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265