

Product Name :
Cylindrical Mould

Product Code :
ASTConcrete-0005



Description :

Cylindrical Mould

Technical Specification :

The cylindrical mould is split vertically into two parts. The mean internal diameter is within ± 0.2 mm and height is within ± 1 mm. The ends are machined to ± 0.05 . The base plate and top plate are machined flat to ± 0.03 mm.

- Mould cylinder, cast iron, 150mm dia x 150mm height.
- Mould cylinder, cast iron, 150mm dia x 300mm height.
- Mould cylinder, cast iron, 160mm dia x 320mm height.
- Mould cylinder, cast iron, 100mm dia x 200mm height.
- Mould cylinder, cast iron, 300mm dia x 600mm height.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265