

**Product Name :**  
Portable Core Cutting Drilling Machine Petrol Driven

**Product Code :**  
ASTROCK-0002



**Description :**

Portable Core Cutting Drilling Machine Petrol Driven

**Technical Specification :**

Specification:

Compact and portable Labtest core drilling machine is designed to cut cores up to 200 mm diameter from concrete, asphalt and similar hard construction materials. The machine comprises a vertical support column which carries the drill head/motor assembly. The motor assembly comprises a 6.5 hp petrol engine. A ball screw mechanism enables close control of the drilling pressure and rapid return when drilling is completed. A water spraying assembly is mounted on the machine. The complete assembly is supplied on a rigid wheel mounted metal base frame with leveling and fixing facility during the operation.

Weight (approx) : 100 Kg

Power : Petrol Engine 6.5 HP

Briggs & Stratton USA



## Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265