

Product Name :
Digital PH Meter

Product Code :
ASTENVIOR-0003



Description :

Digital PH Meter

Technical Specification :

DIGITAL PH METER

SPECIFICATIONS

PH Range 0 To 14 PH

MV Range 0 To +1999 MV

Resolution 0.01 PH, In PH Range 1 Mv In MV Range

Repeatability 0.01 PH 1 Digit 1 MV 1 Digit

Standardization Range DLABn. 2 PH.

Temp Compensation 0 - 1000c

Display 3½ Digits LED

Power 230v 10%, 50 Hz

DIGITAL PH METER

SPECIFICATIONS

PH Range 0 To 14 PH

MV Range 0 To +1999 MV

Resolution 0.01 PH, In PH Range 1 Mv In MV Range

Repeatability 0.01 PH 1 Digit 1 MV 1 Digit

Standardization Range DLABn. 2 PH.

Temp Compensation 0 - 1000c

Display 0 - 1000c Manual And Automatic Compensation With (PT 100)

Power 230v 10%, 50 Hz

DLABCROPROCESSOR PH MODEL

DISPLAY 16 X 2 Alpha Numeric LCD Display.

RANGE pH : 0 To 14.00 mV : 0 +/- 1999.9.

RESOLUTION 0.01 pH, 0.01 mV, Temp. 0.01UC.

ACCURACY $\hat{A}\pm 0.01$ pH, $\hat{A}\pm 1$ mV, Temp. $\hat{A}\pm 1$ UC. (Auto/Manual)

TEMPERATURE COMPENSATION 0 to 100UC

CALIBRATION Auto / Manual.

STORAGE Up to 100 Samples.

POWER 230V $\hat{A}\pm 10\%$, 50 Hz AC

AUTO BUFFER RECOGNIITION Auto 4.00, 7..00, 9.20pH.



Astra Scientific

www.astrascientific.com, Email: info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 Phone: +91-8860605265