

Product Name :
Conductivity Meter

Product Code :
ASTENVIOR-0004



Description :

Conductivity Meter

Technical Specification :

DIGITAL PH METER

DLABCROPROCESSOR CONDUCTIVITY / TDS / TEMP.

DISPLAY 16 X 2 Alpha Numeric LCD Display

RANGE 0.1 μ S to 100mS (6 Decadic Ranges)

RESOLUTION 0.01 μ S.

ACCURACY \pm 1% of F.S. \pm 1 Digit

CALIBRATION Auto / Manual.

TEMPERATURE COMP. 0 to 100 UC.

CELL CONSTANT 0.1 To 5.0 Adjustable.

SENSOR Up to 100 Samples.

POWER 230V $\hat{\pm}$ 50 Hz. AC.

DELUXE DIGITAL CONDUCTIVITY METER

DISPLAY 3.5 Digit LED.

RANGE 5 Ranges, 0 to 1000Mhos / cm.

ACCURACY $\hat{\pm}$ 1% of F.S. $\hat{\pm}$ 1 Digit.

TEMPERATURE COMP. Auto / Manual 0 To 50 deg. C

CELL CONSTANT Adjustable.

SENSOR Conductivity Cell.

POWER 230V $\hat{\pm}$ 50 Hz AC.

DIGITAL CONDUCTIVITY METER

DISPLAY 3.5 Digit LED.

RANGE 5 Ranges, 0 to 1000Mhos / cm.

ACCURACY $\hat{\pm}$ 1% of F.S. $\hat{\pm}$ 1 Digit.

TEMPERATURE COMP. Manual 0 To 50 deg. C

CELL CONSTANT Adjustable.

SENSOR Conductivity Cell.

POWER 230V $\hat{\pm}$ 50 Hz AC.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265