

Product Name :
Slump Cone Test Apparatus

Product Code :
ASTConcrete-0007



Description :

Slump Cone Test Apparatus

Technical Specification :

Ref. Standards IS:7320, BS:1881 , ASTM C 143, AASHTOT119

Slump Cone Test Apparatus for testing consistency and workability of fresh concrete meet ASTM C143, AASHTO T 119 and BS 1881 standards. The slump cone having dimension as (ID-100mm at top ID- 200mm at bottom & length 300mm) is filled with freshly mixed concrete and is tamped with a tamping bar in four layers. The top of the concrete is levelled off with the top of the slump cone, the cone is lifted vertically up and the slump of the sample is immediately measured. Slump Test Apparatus is with testing rod & base plate.

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com