

Product Name :
Laboratory Drum Type Mixer

Product Code :
ASTConcrete-0011



Description :

Laboratory Drum Type Mixer

Technical Specification :

To find out the effects of the different type constituents on the concrete properties requires the preparation of numerous and smaller volume of concrete batches in the laboratory.

Dimensional volume of the mixing pan is 125 liters and the effective mixing capacity is 75 liters. Its ability of mixing small volume of mixtures with high performance through the second engine which rotates the stirring beater in the opposite direction of the mixing drum's rotation direction.

The mixing drum can be tilted to empty the mixture without hassle and it can be removed for easy cleaning on completion of the mixing operation. It is equipped with rubber wheels which provides high mobility.

All parts of the mixer are galvanized and painted with non-corrosive paint. The protection class of machine is IP55.

Drum volume: 125 litres

Yield: 75 litres of concrete

Power supply: 230V 1ph 50 Hz - 0,3 HP

Dimensions: 720x1320x1280 mm

Weight: 60 Kg



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265