

Product Name :
V-Funnel Test Apparatus

Product Code :
ASTConcrete-0024



Description :

V-Funnel Test Apparatus

Technical Specification :

General Properties:

V Funnel apparatus is used to evaluate the flow time of freshly mixed self-compacting concrete. The test is not suitable when the maximum size of the aggregate exceeds 22.4 mm.

The test set consists of a stainless steel funnel placed vertically on a supporting stand. The discharge orifice is equipped with a lid, which can be momentarily opened.

The L Shape Box Apparatus is supplied complete with;

- Filling Hopper
- Base

Dimensions 525x300x1040 mm



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265