


Product Name :
Pocket Concrete Penetrometer

Product Code :
ASTConcrete-0039



Description :

Pocket Concrete Penetrometer

Technical Specification :

Specification:

Consists of a needle having face area $3/10$ sq. cm. and graduated at a distance of 25 cm. The needles point is an integral part of barrel which houses a calibrated spring. The spring is confined in a sleeve. The resistance offered by the concrete mortar is shown on the direct reading scale with a marker ring which holds its position when released. Scale range is 0-50 kg/cm² when the penetration resistance reaches a value of 35kg/cm² the concrete is assumed initially set. Supplied complete in carrying case.

Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265

