

Product Name :
Bleeding Of Freshly Mixed Concrete

Product Code :
ASTConcrete-0040



Description :

Bleeding Of Freshly Mixed Concrete

Technical Specification :

Specification: This is used for determination of the relative quantity of mixing water that will bleed from a sample of freshly mixed concrete. Stainless steel cylindrical container for bleeding test to EN 480-4, approx. 14 liters capacity.

Internal dia 255 +/- 5mm x 280 +/- 5mm inside height. Complete with stainless steel cover. Flat rounded steel float with a diameter of (100 ± 10) mm.



Astra Scientific

www.astrascientific.com, Email: info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 Phone: +91-8860605265

