

Product Name :
Unconfined Compression Testing Apparatus

Product Code :
ASTSOIL-0012



Description :

Unconfined Compression Testing Apparatus

Technical Specification :

Specification:

Comprises a screw operated load frame, cap.5000kg with a gear box and motor drive giving 1.25, 1.5, 2.5mm/min. Rates of strain, a pair of cone seating, adaptor for proving ring, and strain dial gauge bracket. Supplied with one pair of male/female coning tools for 38mm dia samples but without ring and dial gauge. Suitable for operation on 230 V, Single Phase A.C. Supply.

Accessories:

Coning tools in pairs (Male & Female) for samples having diameter 38, 50, 75 & 100mm



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265