

Product Name :
RECIPROCATING PUMP DEMONSTRATION UNIT, Computer
Interface

Product Code :
FM177



Description :

A bench top unit with a reciprocating pump, 0.37kW motor with reduction gear, storage tank measuring tank, pressure relief valve, pulse chambers, sensors, and software. To be used with HP 100, HP 107 and computer or HP 104 and computer.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265