


Product Name :
Proctor penetrometer

Product Code :
ASTSOIL-0029



Description :

Proctor penetrometer

Technical Specification :

Specification:

Used for establishing the moisture content-penetration resistance relationship of fine-grained soils. It consists of a special spring dynamometer with pressure indicating scale on the stem of the handle. A sliding ring on the stem indicate the maximum pressure obtained in the test. Supplied in a wooden carrying case.

- Load scale: 0 to 55 kg , 1 kg subdivisions with max load indicator
- Diameter of interchangeable needles: 28.55, 24.79, 20.22, 16.54, 12.83, 9.07, 6.40, 5.23 and 4.52 mm



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265