


Product Name :
Hydraulic Shrinkage Mould 25x25x285 mm

Product Code :
ASTCEMENT-0006



Description :

Hydraulic Shrinkage Mould 25x25x285 mm

Technical Specification :

General Properties:

Two gang prism mould 25x25x285 mm complete with contact points. Gauge length 250mm, total length 285 mm.

Weight approx.: 6 kg. This mould is the same used for the Expansion of Portland cement with the High pressure autoclave, conforming to ASTM C151 and ASTM C490.

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265



Astra Scientific