

Product Name :

Cement - Autoclave, high pressure: Autoclave, high pressure

Product Code :

ASTCEMENT-0007

**Description :**

Cement - Autoclave, high pressure: Autoclave, high pressure

Technical Specification :**General Properties:**

The autoclave consists of a high-pressure steam vessel with internal dimensions 154 mm dia. x 430 mm high to accept a rack for holding 10 specimens obtained with the moulds (see accessories). Complete with pressure gauge, pressure regulator, temperature regulator, control switches, safety valve and specimen rack.

Certified conforming to ISPELS procedure.

- Heater: 2600 w
- Overall dimensions: 450x475x1080 mm



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265