


**Product Name :**  
Two gang Prism Mould 25x25x285 mm

**Product Code :**  
ASTCEMENT-0008



**Description :**

Two gang Prism Mould 25x25x285 mm

**Technical Specification :**

**General Properties:**

This mould is the same used for the Expansion of Portland cement with the High pressure autoclave, conforming to ASTM C151 and ASTM C490. Two gang prism mould 25x25x285 mm complete with contact points.



**Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265