


Product Name :
Two gang Prism Mould 25x25x285 mm

Product Code :
ASTCEMENT-0008



Description :

Two gang Prism Mould 25x25x285 mm

Technical Specification :

General Properties:

This mould is the same used for the Expansion of Portland cement with the High pressure autoclave, conforming to ASTM C151 and ASTM C490. Two gang prism mould 25x25x285 mm complete with contact points.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265