


Product Name :
Prism Mould 25x25x280 mm to UNI

Product Code :
ASTCEMENT-0009



Description :

Prism Mould 25x25x280 mm to UNI

Technical Specification :

General Properties:

Used for accelerated expansion test on mortar prisms to UNI 11504. Three gang prism mould 25x25x280 mm complete with contact points. Gauge length 250mm, total length 280 mm.

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com