

Product Name :
Prism Mould 75x75x285 mm to UNI

Product Code :
ASTCEMENT-0011



Description :

Prism Mould 75x75x285 mm to UNI

Technical Specification :

General Properties:

Used for expansion test on concrete prisms to UNI 11604. One gang prism mould 75x75x285 mm complete with contact points. Gauge length 250mm, total length 285 mm.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265