

Product Name :
Steel Measure, 400 ml

Product Code :
ASTCEMENT-0012



Description :

Steel Measure, 400 ml

Technical Specification :

General description:

Used as a density method for the determination of the air content of freshly mixed mortars. Made of steel, 88.1 mm height x 76.2 mm internal diameter, it is calibrated to contain 400 ± 1 ml of water at 23°C.



Astra Scientific

www.astrascientific.com, Email: info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265