


Product Name :
Pat Rest, Brass Ring Mould

Product Code :
ASTCEMENT-0013



Description :

Pat Rest, Brass Ring Mould

Technical Specification :

General Properties:

Used for the determination of the soundness of hydrated lime and gypsum plasters. The mould is 100 mm dia. x 5 mm deep with an internal taper of 5°. It is supplied complete with glass plate. Three are required for each test.

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com