

Product Name :
Water Retention Apparatus

Product Code :
ASTCEMENT-0017



Description :

Water Retention Apparatus

Technical Specification :

General Properties:

Used for determining the water retention value of cement and lime.

Two models available:

- one fitted with aspirator pump and analog vacuum gauge with regulator

Water retention apparatus including an aspirator pump, vacuum regulator, vacuum gauge, three way stopcock, flask, rubber gasket, glass funnel, perforated brass dish and stand.

Weight: 8 kg (approx.)

- the other with a portable vacuum pump and digital vacuum gauge with regulator

Water retention apparatus including vacuum pump with ultimate vacuum 0.1mbar, digital vacuum regulator complete with vacuum gauge 1mbar resolution, three way stopcock, flask, rubber gasket, glass funnel, perforated brass dish and stand. 230 V/50-60Hz/1ph



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265